

Data recovery services

Mosaiq Lantis project since 2018
BeamView Images

Initial situation:

Information about a patient is encoded and distributed in different directories in Mosaiq. That can be, for example, DRRs or CT data. What makes it more difficult is that the files have to be prepared differently depending on the manufacturer and software version. BeamView images also need to be adjusted.

A simple transfer to a DICOM archive is not possible.

Our solution:

Tiger-Software has developed various tools to prepare old data from the Mosaiq database sets into a valid DICOM format.

End goal:

Storage in a DICOM archive

Customer raw files from Mosaiq (Lantis) are to be decoded and modified into RT-IMAGE or Secondary Capture image according to the standard depending on the recognized system. Additional conversion of patient photos into DICOM Secondary Capture Images, JPG or PDF files is also possible.

Alternative: Local archive











It is possible to send the files to other systems via ‘DICOM send’ by entering the patient ID or searching for the name.

One option is to display the DICOM images with our software TigClick.


What do we need to do

There are different modality's in different Mosaiq directory's (BDS, Image, Image1)
CT, RTDOSE, RTPLAN, RTSTRUCT, RTIMAGE, OT, SC, Patient-Photo

For each manufacturer and modality we have to check and adapt the files.

 BEAMVIEW COHERENCE	 Pinnacle3
 COHERENCE Therapist Workspace	 SIEMENS
 IMPAC Medical Systems, Inc	 Siemens AG
 iPlan RT Dose	 syngo Suite for Oncology
 MOSAIQ	 XiO

The result for the user



The processed data records are stored by Tiger-Software in a patient ID-directory.

Based on the fact that discrepancies can occur in the master data, Tiger-Software **checks all files** of each patient with the same ID


Possible errors:

Transposed numbers in the date of birth,
differences in gender,
spelling error in the name

Rework is possible.

Transfer to the archive:

Modified files from directories can be automatically exported via 'DICOM send'.



This layout is available in German or English. Other languages on request.

Contact via

www.tiger-software.de or info@tiger-software.de